

Introduction to ATAC-seq

Shamith Samarajiwa


University of Cambridge

CRUK Summer School in Bioinformatics

July 2018


Higher-order

Primary-order




Techniques	Validation				MNase-seq	DNase-seq	FAIRE-seq	ATAC-seq
	4C	5C	Hi-C	ChIP, Imaging etc.				
	1-vs-Many	Many-vs-Many	All-vs-All		Inferred closed regions	Open chromatin	Open chromatin	Open/Closed chromatin

Procedure	Higher-order		Primary-order	
	Experimental	Computational	Experimental	Computational
	1. Crosslinking	1. Alignment	1. Crosslinking*	1. Size-selection
	2. Fragmentation	2. Filtering	2. Fragmentation	2. Alignment
	3. Ligation	3. Binning	3. Size-selection	3. Peak calling
	4. Detection	4. Normalization	4. Sequencing	4. Normalization
		5. Identifying interactions		5. Visualization
		6. Visualization		




Assay for Transposase Accessible Chromatin






ATAC-seq

- Measure **chromatin accessibility** and **nucleosome free regions**
- Can also detect nucleosome packing, positioning and TF footprints
- Does not require sonication and phenol chloroform extractions, antibodies or sensitive enzymatic digestions that can introduce potential bias
- A hyperactive Tn5 transposase is used to fragment DNA and integrate into active regulatory regions
- During ATAC-seq, 500–50,000 unfixed nuclei are tagged *in vitro* with sequencing adapters by Tn5 transposase




ATAC-seq

- Two-step protocol
 - Insertion of Tn5 transposase with adaptors
 - PCR amplification
- Needs ~500-50,000 cells
- Paired-end reads produce information about nucleosome positioning
- Insert size distribution of fragments has a periodicity of ~200 bp, suggesting that fragments are protected by multiples of nucleosomes
- Different fragmentation patterns can be associated with different functional states (eg. TSSs are more accessible than promoter flanking or transcribed regions)



[1] JD Buenrostro et al, *Nature Methods*, 2013.

Workflow of ATAC-seq data processing



Working with ATAC-seq

- Remove mitochondrial reads
 - A large fraction of ATAC-seq reads map to the mitochondrial genome (up to 40-60%)
- Remove blacklisted regions before peak calling
 - hg19 or hg38 blacklisted regions can be obtained from ENCODE
- The signal (open chromatin where the transposase was inserted) is a mixture of feature types:
 - Short fragments - signal from nucleosome free regions (**NFRs**) and open regions around DNA bound transcription factors. TSS rich.
 - Longer fragments - open regions from around nucleosomes. Includes +1 and -1 nucleosome positions

Differences from ChIP-seq data processing

- Peak calling: Use fragment size for smoothing when calling peaks with MACS2
- Alternatively, to identify open regions following MACS parameters are used:
MACS2 callpeak -t bamfile --nomodel --shift -100 --extsize 200 --format BAM -g hg38

For nucleosome occupancy shift and extension can centre the signal on nucleosomes (147 bp DNA is wrapped in a nucleosome)

MACS2 callpeak -t bamfile --nomodel --shift 37 --extsize 73 --format BAM -g hg38

Optionally, for paired end data

MACS2 callpeak -t bamfile --format BAMPE -g hg38

- --nomodel: don't build shifting model
- --shift: when this value is negative, ends will be moved toward 3'→5' direction
- --extsize: extend reads in 5'→3' direction to fix-sized fragments

Normalization and Differential Accessibility of ATAC-seq

- Normalisation across samples might be needed
 - Efficiency of the ATAC-seq protocol in assaying open regions is affected by how many transposons get into the nuclei
 - One normalisation solution is to use signal from ‘essential or housekeeping genes’. see: *Sarah K. Denny et al, Cell, 2016.*
- **THOR** is an HMM-based approach to detect and analyze differential peaks in two sets of ChIP-seq data from distinct biological conditions with replicates
 - Normalizes bam files (given a bed file of housekeeping genes)
 - Can compare two ATAC-seq datasets to perform differential openness analysis

Cutting site and TF Footprinting analysis

- ATAC-seq produces shorter fragments around smaller protected regions such as TF bound regions.
- Cut sites pileup signal is a good guide to TF binding sites
- Shifting reads +4/-5 nt depending on strand, should adjust for expected shift from Tn5 insertion.
- TF footprinting of these regions enables the detection of motifs for bound TFs
 - Centipede
 - msCentipede
 - PIQ
 - Mocap